

SIT support

**Innovations for
your dental practice**

BODY support



The solution for treatment with long anaesthetic time

- Optimal position for treatment via perfect support for the body
- 7 cm high visco-elastic foam core with memory effect
- Safe anchorage on the treatment unit via anti-slip fittings at the bottom
- Comfortable to lie on due to high quality artificial leather cover with removable foot rest

Perfect positioning of the “small” patients

- Avoids early signs of tiredness via comfortable cushion elements with unique lower leg supports
- Safe anchorage on every treatment unit via flexible adjustment and antislip fittings at the bottom
- Hard-wearing, dirt-resistant surface made of high quality, permanently elastic and anti allergic artificial leather cover

KID support



HEAD & NECK support



Provides ergonomically correct and relaxing support

- Perfect anchorage for the head and optimal treatment via visco-elastic foam core and a “memory-effect”
- No tension in the neck muscles via pressure relief and ergonomic design
- Comfortable to lie on due to high quality artificial leather cover

M&C PLUS support

Suitable for treatment units with concertina joint on the adjustable headrest

- Perfect anchorage for the head via visco-elastic foam core with a “memory-effect” reduces spontaneous movement and allows optimal treatment
- Ergonomic design ensures a relaxed head position which promotes blood flow and prevents tension in the neck muscles
- Cushioned “plus” extension protects the neck region from pressure points
- Soft pleasant breathable cover gives your patients the highest comfort
- Particularly easy to remove and clean with standard disinfection agents



- Free accessible operating elements